

## Cold light source

### Compact dual 2 W LED

#### HIGHLIGHTS

- Economical and compact light source
- For Materials Science en stereo microscopy
- 2x 1 W LED
- With two flexible gooseneck arms
- 2 years warranty



● LE.5208

#### TECHNICAL SPECIFICATIONS

- Economical dual LED intensity adjustable light source
- Equipped with two 44 cm flexible gooseneck arms, each equipped with a 1 W LED
- The light source version does not require bulb changing like halogen light sources and has a longer life span and low energy consumption
- The lightsource produces white light and there is no color temperature shifting when dimming
- Ideal for life and materials sciences, various industrial applications and for stereo microscopy

#### MODEL

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Power</b>             | 2x 1 W LED (equivalent to 2x 20 W halogen) |
| <b>Luminous flux</b>     | 18,000 lux max. (2 LED's)                  |
| <b>Color temperature</b> | 6,500 K                                    |
| <b>Operating voltage</b> | 100-240 V AC / 12 V DC (50/60 Hz)          |
| <b>Dimensions</b>        | (h) 70 x (w) 120 x (d) 170 mm              |
| <b>Weight</b>            | 1.6 kg                                     |
| <b>Product Number</b>    | LE.5208                                    |