

Matrix for KITOZOOM

Adapter 0,5x		attachement lense 0,25X	attachement lense 0,5x	attachement lense 0,75x	without 1,0x	attachement lense 1,5x	attachement lense 2x
Zoom position 0,7	field of view mm(horizontal)	56	28.1	18.4	12.7	9	5.2
	magnification (24" monitor)	8 x	15 x	23 x	34 x	48 x	83 x
Zoom position 1	field of view mm (horizontal)	37	18.5	12.1	9	6	4.4
	magnification (24" monitor)	12 x	23 x	36 x	48 x	72 x	98 x
Zoom position 2	field of view mm(horizontal)	18.5	9.2	6	4.4	2.95	2.2
	magnification (24" monitor)	23 x	47 x	72 x	98 x	146 x	195 x
Zoom position 3	field of view mm (horizontal)	12	6.1	4	2.9	1.95	1.45
	magnification (24" monitor)	36 x	70 x	107 x	148 x	220 x	297 x
Zoom position 4	field of view mm(horizontal)	9	4.6	3	2.2	1.45	1.05
	magnification (24" monitor)	48 x	93 x	143 x	195 x	297 x	410 x
Zoom position 4,5	field of view mm (horizontal)	8	4	2.65	1.95	1.3	0.95
	magnification (24" monitor)	54 x	107 x	162 x	220 x	331 x	453 x
working distance		310 mm	178 mm	116 mm	89 mm	47 mm	27 mm
system length		243	235	235	230	236	239

Note: Standard configuration KITOZOOMs include only 0.5x attachment lense, others are optional